



BS-DP Biliary Stents Double Pigtail



The BS-DP Biliary Stents Double Pigtail are used for drain obstructed biliary ducts.

- Radiopaque polyurethane tubing material
- Taper tip for easy insertion
- Smooth rounding posterior
- Various sizes available



Order No.	Fr	Length (cm)	Qty. / Box
BS-DP5/4	5	4	10
BS-DP5/7	5	7	10
BS-DP5/10	5	10	10
BS-DP7/4	7	4	10
BS-DP7/7	7	7	10
BS-DP7/10	7	10	10
BS-DP10/4	10	4	10
BS-DP10/7	10	7	10
BS-DP10/10	10	10	10

Order No.	Fr	Length (cm)	Qty. / Box
BS-DP 5/4-K	5	4	1
BS-DP 5/7-K	5	7	1
BS-DP 5/10-K	5	10	1
BS-DP 7/4-K	7	4	1
BS-DP 7/7-K	7	7	1
BS-DP 7/10-K	7	10	1
BS-DP 10/4-K	10	4	1
BS-DP 10/7-K	10	7	1
BS-DP 10/10-K	10	10	1

*All products may not be available in all countries. Please contact our sales department for details.

Kit includes :

- Biliary Stents
- Majestick - Biliary Stent Introducer
- 0.035" PTFE Guidewire 400cm