



BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip




The BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic band on the inner catheter to provide visibility under ultrasound
- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe and outer sheath
- Bulb shaped angled tip of the outer sheath
- Silicone depth positioner for precise placement
- The catheter is curved to facilitate insertion



Order No.	Description	Qty. / Box
BT-ETC-U-S	BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip	10
STY-BU	Stylet for BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip	10

*All products may not be available in all countries. Please contact our sales department for details.



Inner catheter tip
echogenic metal band