



BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip





The BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe
- Bulb shaped angled tip of the outer sheath
- Silicone depth positioner for precise placement
- The bulb-shaped angled outer sheath eases passage through the cervix
- Echogenic band on the inner catheter to improve the visibility under ultrasound



Order No.	Description	Qty. / Box
BT-ETC-U-S-P	BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip	10
STY-BU	Stylet for BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip	10

*All products may not be available in all countries. Please contact our sales department for details.



Inner catheter tip
echogenic metal band