



CURVED Intra Uterine Insemination Catheter



The CURVED Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Curve shaped catheter
- Open tip as well as closed tip (with two lateral holes) catheters available



Order No.	Description	Qty. / Box
I-IUI-C	CURVED Intra Uterine Insemination Catheter	25
I-IUI-C-O	CURVED Intra Uterine Insemination Catheter (Open Tip)	25

*All products may not be available in all countries. Please contact our sales department for details.



Close Tip

Open Tip