



DUALFLOW Double Lumen Ovum Pickup Needle



The DUALFLOW OVUMPICK Double Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- The outer lumen allows flushing of follicles



Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
DUF-DL17/35	DUALFLOW Double Lumen Ovum Pickup Needle	17	35	10

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip