



FALCON Single Lumen Ovum Pickup Needle




The FALCON Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound - guided transvaginal retrieval of oocytes from ovarian follicles.

- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line



Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
F-OVM-SL17/35	FALCON Single Lumen Ovum Pickup Needle	17	35	10

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip