


FLEXI Intra Uterine Insemination Catheter



The FLEXI Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- Shapeable outer sheath for easy insertion
- Centimetre marking on a catheter for precise placement
- Stylet available separately
- Open tip as well as closed tip (with two lateral holes) catheters available

Centimetre marking provided



Soft uterine catheter



Shapeable
outer sheath


Firm
support



Stylet

Order No.	Description	Qty. / Box
F-IUI	FLEXI Intra Uterine Insemination Catheter	25
F-IUI-O	FLEXI Intra Uterine Insemination Catheter (Open Tip)	25
STY18	Stylet for FLEXI IUI Intra Uterine Insemination Catheter	25

*All products may not be available in all countries. Please contact our sales department for details.



Close Tip



Open Tip