


## FlexiGuide Standard - J Tip Guidewire



The FlexiGuide Standard - J Tip Guidewire is used for percutaneous entry into vessel using the seldinger technique or for endoscopic access.

- J Tip fixed core
- PTFE coating allows smooth catheter insertion and stent tracking



Order No.	Tip	Diameter	Length (cm)	Qty. / Box
JFT150-25-3	3mm, J-Tip Fixed Core, PTFE Coated	0.025"	150	5
JFT150-32-3	3mm, J-Tip Fixed Core, PTFE Coated	0.032"	150	5
JFT150-35-3	3mm, J-Tip Fixed Core, PTFE Coated	0.035"	150	5
JFT150-38-3	3mm, J-Tip Fixed Core, PTFE Coated	0.038"	150	5

\*All products may not be available in all countries. Please contact our sales department for details.