




HYSTERO CATH Hysterosalpingography Catheter



The HYSTERO CATH Hysterosalpingography Catheter is used for injection of contrast media into the uterine cavity and fallopian tube(s) to evaluate the patency of the fallopian tubes.

- Strong low profile 7 Fr balloon
- Closed tip with side hole and open tip HSG catheter available
- Stopcock designed for locking the balloon channel and regulating the flow of contrast media.
- Stylet provided for added rigidity during the introduction
- 3 ml syringe included



Order No.	Description	Qty. / Box
HC-HSG7	HYSTERO CATH Hysterosalpingography Catheter	5
HC-HSG-O7	HYSTERO CATH Hysterosalpingography Catheter Open Taper Tip	5

*All products may not be available in all countries. Please contact our sales department for details.

Caution:

This product contains natural rubber latex which may cause allergic reaction in some individuals.