



INTEGRA Intra Uterine Insemination Catheter



The INTEGRA Intra Uterine Insemination Catheter (Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- Placement becomes easier with the distal tip's soft, flexible curve.
- The stainless steel inner cannula is malleable.
- Luer lock connector.
- Silicone stopper for adjustable positioning.



Order No.	Description	Qty. / Box
IG-IUI	INTEGRA Intra Uterine Insemination Catheter	25

*All products may not be available in all countries. Please contact our sales department for details.



Side Hole