


OPICK Single Lumen Ovum Pickup Needle with Knurling



The OPICK Single Lumen Ovum Pickup Needle with Knurling is used for laparoscopic or ultrasound - guided transvaginal retrieval of oocytes from ovarian follicles.

- Echogenic tip with knurling
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line

Tri-facet geometry
with a echogenic tip

Ergonomic handle

Aspiration line


Silicone bung
for tight fit

Luer lock connector

Vacuum line

Order No.	Description	Needle Gauge	Needle Length (cm)	Qty. / Box
OPN-SL17/35-K	OPICK Single Lumen Ovum Pickup Needle with Knurling	17	35	10

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip
with knurling