



OVUMPICK Double Lumen Ovum Pickup Needle



The OVUMPICK Double Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- The outer lumen allows flushing of follicles



| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|-------------|--|--------------|--------------------|------------|
| OVM-DL16/35 | OVUMPICK Double Lumen Ovum Pickup Needle | 16 | 35 | 10 |
| OVM-DL17/35 | OVUMPICK Double Lumen Ovum Pickup Needle | 17 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip