



PRO ECHO TRANS Embryo Transfer Catheter




The PRO ECHO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic inner catheter facilitates enhanced visibility under ultrasound
- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Silicone depth positioner for precise placement
- Orientation arrow on hub
- Angled outer sheath provides easy passage through the cervical canal
- Available in compatible stylet or soft obturator



Order No.	Description	Qty. / Box
PET-ETC	PRO ECHO TRANS Embryo Transfer Catheter	10
PET-ETC-SB	PRO ECHO TRANS Embryo Transfer Catheter with Soft Obturator	10
PET-ETC-STY	PRO ECHO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter	10
STY-PE	Stylet for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.



**Full echogenic
inner catheter**