



PRO EMBRYO STAR Embryo Transfer Catheter



The PRO EMBRYO STAR Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic band on the inner catheter to provide visibility under ultrasound
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Silicone depth positioner for precise placement
- Orientation arrow on the hub
- Angled outer sheath provides easy passage through the cervical canal
- Available in compatible stylet or soft obturator




Soft Obturator for PRO EMBRYO STAR Embryo Transfer Catheter



Stylet for PRO EMBRYO STAR Embryo Transfer Catheter

Order No.	Description	Qty./Box
PEB-ETC-S	PRO EMBRYO STAR Embryo Transfer Catheter	10
PEB-ETC-S-SB	PRO EMBRYO STAR Embryo Transfer Catheter with Soft Obturator	10
PEB-ETC-S-STY	PRO EMBRYO STAR Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter	10
STY-PE	Stylet for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter	10



Echogenic metal band