



PRO EMBRYO TRANS Embryo Transfer Catheter



The PRO EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Silicone depth positioner for precise placement
- Orientation arrow on the hub
- Angled outer sheath provides easy passage through the cervical canal
- Available in compatible stylet or soft obturator



Soft Obturator for PRO EMBRYO TRANS Embryo Transfer Catheter



Stylet for PRO EMBRYO TRANS Embryo Transfer Catheter

Order No.	Description	Qty./Box
PEB-ETC	PRO EMBRYO TRANS Embryo Transfer Catheter	10
PEB-ETC-SB	PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator	10
PEB-ETC-STY	PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter	10
STY-PE	Stylet for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.



Supported inner catheter
with smooth tip