



PRO PASS ECHO TRANS Embryo Transfer Catheter



The PRO PASS ECHO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Orientation arrow on the hub
- Available in compatible stylet or soft obturator



Order No.	Description	Qty. / Box
PPE-ETC	PRO PASS ECHO TRANS Embryo Transfer Catheter	10
PPE-ETC-SB	PRO PASS ECHO TRANS Embryo Transfer Catheter with Soft Obturator	10
PPE-ETC-STY	PRO PASS ECHO TRANS Embryo Transfer Catheter with Stylet	10
SB-PE	Soft Obturator for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter	10
STY-PE	Stylet for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.

