


SURELOCK Vitrification Straws



The SURELOCK Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- Concave tip for exact and secure placement of specimen under microscope
- Secure air tight coupling
- Semi-transparent body and cap
- The arrow on top surface ensuring orientation during use
- Easy visibility under microscope
- 1 - cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of 1-cell stage of embryos
- Available in 5 colours



Colour	Pink	Yellow	Green	Blue	Violet	Qty. / Box
Order No.	SUL-P	SUL-Y	SUL-G	SUL-B	SUL-V	10

*All products may not be available in all countries. Please contact our sales department for details.

