


SURESTRIP Vitrification Straws - Weighted



The SURESTRIP Vitrification Straws-Weighted is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- Weighted balls prevents it from floating in nitrogen
- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- 1 - cell MEA $\geq 80\%$ expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours

Ring marking for easy visualization during insertion

Transparent cap

Weighted balls



Marking on tip for easy visibility under microscope



Flat surface wide body for writing

Strip width ensuring minimum placement of solution

Colour	White	Red	Ivory	Blue	Yellow	Qty. / Box
Order No.	SUS-W-WT	SUS-R-WT	SUS-I-WT	SUS-B-WT	SUS-Y-WT	10

*All products may not be available in all countries. Please contact our sales department for details.



Weighted balls