


THE SAMPLER Endometrial Biopsy Catheter



THE SAMPLER Endometrial Biopsy Catheter is used for performing histological tissue biopsy and extracting samples of uterine menstrual content. Used for obtaining endometrial tissue.

- Centimetre marking provided on the catheter for precise placement
- Smooth rounded tip

Smooth rounded tip



Centimetre marking



Order No.	Description	Qty. / Box
SEBC	THE SAMPLER Endometrial Biopsy Catheter	25

*All products may not be available in all countries. Please contact our sales department for details.