



TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet



The TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Malleable stylet to aid catheter insertion
- Open tip as well as closed tip (with two lateral holes) catheters available



Order No.	Description	Qty. / Box
T-IUI	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet	25
T-IUI-O	TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip)	25

*All products may not be available in all countries. Please contact our sales department for details.



Close Tip



Open Tip