



Trial Transfer Catheter for Pro Pass



The Trial Transfer Catheter for Pro Pass is a closed-ended catheter that is only to be used to determine the passage through the cervix prior to embryo transfer accessible to Allwin's Embryo Transfer Catheter.

- Allows for confident positioning of the outer sheath before introducing the embryo loaded catheter
- Same soft material and smooth tip as catheter to minimize patient trauma
- Centimetre depth marking for confident and precise catheter placement
- Limulus Amebocyte Lysate (LAL) tested to ensure the product is non-pyrogenic



Order No.	Description	Qty. / Box
TT-PP-ETC	Trial Transfer Catheter for Pro Pass Embryo Transfer Catheter	10
TT-PPS-ETC	Trial Transfer Catheter for Pro Pass Star Embryo Transfer Catheter	10
TT-PPE-ETC	Trial Transfer Catheter for Pro Pass Echo Embryo Transfer Catheter	10
SB-PE	Soft Obturator for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter	10
STY-PE	Stylet for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter	10

*All products may not be available in all countries. Please contact our sales department for details.

