



WOMEN'S HEALTH

| | |
|--|-------|
| INTRODUCTION | 2 |
| OVUM PICKUP NEEDLES | 3-11 |
| ACE Single Lumen Ovum Pickup Needle..... | 3 |
| ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line..... | 4 |
| EAGLE Single Lumen Ovum Pickup Needle..... | 5 |
| OPICK Single Lumen Ovum Pickup Needle with Knurling..... | 6 |
| FALCON Single Lumen Ovum Pickup Needle..... | 7 |
| ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush..... | 8 |
| ACE-M Manual Single Lumen Ovum Pickup Needle..... | 9 |
| DUALFLOW Double Lumen Ovum Pickup Needle..... | 10 |
| OVUMPICK Double Lumen Ovum Pickup Needle..... | 11 |
| EMBRYO TRANSFER CATHETERS | 12-25 |
| EMBRYO TRANS Embryo Transfer Catheter..... | 12 |
| ECHO TRANS Embryo Transfer Catheter..... | 13 |
| Trial Transfer Catheter..... | 14 |
| BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip..... | 15 |
| BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip..... | 16 |
| BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip..... | 17 |
| BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip..... | 18 |
| PRO EMBRYO TRANS Embryo Transfer Catheter..... | 19 |
| PRO EMBRYO STAR Embryo Transfer Catheter..... | 20 |
| PRO ECHO TRANS Embryo Transfer Catheter..... | 21 |
| PRO PASS Embryo Transfer Catheter..... | 22 |
| PRO PASS STAR Embryo Transfer Catheter..... | 23 |
| PRO PASS ECHO TRANS Embryo Transfer Catheter..... | 24 |
| Trial Transfer Catheter for Pro Pass..... | 25 |
| INTRA UTERINE INSEMINATION CATHETERS | 26-31 |
| INSEMINA Intra Uterine Insemination Catheter..... | 26 |
| CURVED Intra Uterine Insemination Catheter..... | 27 |
| TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet..... | 28 |
| FLEXI Intra Uterine Insemination Catheter..... | 29 |
| INTEGRA Intra Uterine Insemination Catheter..... | 30 |
| STUMPY Intra Uterine Insemination Cannula..... | 31 |
| FALLOPIAN TUBE RECANALIZATION CATHETER | 32 |
| SALPINX CATH Fallopian Tube Recanalization Catheter..... | 32 |
| HYSTEROSALPINGOGRAPHY CATHETER | 33-34 |
| HYSTERO CATH Hysterosalpingography Catheter..... | 33 |
| HYSTERO CATH SILICONE Hysterosalpingography Catheter..... | 34 |
| ENDOMETRIAL BIOPSY SAMPLING SYSTEM | 35 |
| THE SAMPLER Endometrial Biopsy Catheter..... | 35 |
| PRO-VAC Endometrial Biopsy Sampler..... | 35 |
| VITRIFICATION STRAWS | 36-38 |
| SURESTRIP Vitrifcation Straws..... | 36 |
| SURESTRIP Vitrifcation Straws - Weighted..... | 37 |
| SURESTRIP Vitrifcation Straws - Weighted Plug..... | 38 |
| SURELOCK Vitrifcation Straws..... | 39 |
| D-TIP Denudation Pipettes..... | 40 |
| AspiPump 90D..... | 41 |
| OTHER PRODUCTS | 42 |
| ASPI Aspiration Needle..... | 42 |

Product Warranty

Allwin warrants that at the time of manufacturing, the products are prepared and tested in accordance with the highest standards. Due to biological differences in individuals, no product is 100% effective under all circumstances for all patients. Further, because we have no control of the condition under which the product is used, diagnosis of the patient, the method of use or administration & handling and storage of the product after it leaves our possession, we do not warrant any consequence which may arise about the final use of the device. The foregoing warranty is exclusive and instead of all other written, oral or implied warranties. No representative of **Allwin** may change any of the foregoing and the buyer accepts the product subject to all terms hereof.

Design Changes

Continuous improvements is the chosen path at **Allwin** and we reserve the right to change or discontinue any product without notice.

Return Goods Policy

Allwin will accept returned goods only after issuing a Return Goods Authorisation (RGA). The return request must provide a reason for return, invoice number and order number. Credit will be issued once the request is evaluated as per company rules. We will accept a return of the stock item within 30 days of the original invoice date. We will not accept the return of a nonstock, opened and damaged products. Returned stock items may be subject to a restocking fee.


Product testing

Allwin products are endotoxin tested with a limulus amebocyte lysate (LAL) assay unless otherwise noted. Our Ovum Pickup needles and Embryo Transfer catheters are batch tested for toxicity with a one-cell mouse embryo assay (MEA) and our Intra Uterine Insemination catheters are Human Sperm Survival Assay (HSSA) tested.

Storage of Products

Allwin recommends that products be stored in a cool and dry place. Plastic material used may degrade following extended periods of exposure to light and temperature.

ACE Single Lumen Ovum Pickup Needle



The ACE Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line

Tri-facet geometry
with a echogenic tip

Ergonomic handle

Aspiration line

Silicone bung
for tight fit

Luer lock connector

Vacuum line


| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|---------------|-------------------------------------|--------------|--------------------|------------|
| A-OVM-SL16/35 | ACE Single Lumen Ovum Pickup Needle | 16 | 35 | 10 |
| A-OVM-SL17/35 | ACE Single Lumen Ovum Pickup Needle | 17 | 35 | 10 |
| A-OVM-SL18/35 | ACE Single Lumen Ovum Pickup Needle | 18 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Echogenic tip

ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line




The ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Aspiration and flushing line is made possible through a single channel
- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line
- Multiple operation functionality by using stopcock valve




| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|------------------|---|--------------|--------------------|------------|
| A-OVM-SL16/35-FL | ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line | 16 | 35 | 10 |
| A-OVM-SL17/35-FL | ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line | 17 | 35 | 10 |
| A-OVM-SL18/35-FL | ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line | 18 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip

EAGLE Single Lumen Ovum Pickup Needle



The EAGLE Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Sharp tip for easier penetration
- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line

Tri-facet geometry
with a echogenic tip

Ergonomic handle

Aspiration line

Silicone bung
for tight fit

Luer lock connector

Vacuum line


| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|---------------|---------------------------------------|--------------|--------------------|------------|
| E-OVM-SL17/35 | EAGLE Single Lumen Ovum Pickup Needle | 17 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Echogenic tip


OPICK Single Lumen Ovum Pickup Needle with Knurling



The OPICK Single Lumen Ovum Pickup Needle with Knurling is used for laparoscopic or ultrasound - guided transvaginal retrieval of oocytes from ovarian follicles.

- Echogenic tip with knurling
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line

Tri-facet geometry
with a echogenic tip




| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|---------------|---|--------------|--------------------|------------|
| OPN-SL17/35-K | OPICK Single Lumen Ovum Pickup Needle with Knurling | 17 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Echogenic tip
with knurling

FALCON Single Lumen Ovum Pickup Needle




The FALCON Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound - guided transvaginal retrieval of oocytes from ovarian follicles.

- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Ergonomic handle
- 100 cm aspiration line




| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|---------------|--|--------------|--------------------|------------|
| F-OVM-SL17/35 | FALCON Single Lumen Ovum Pickup Needle | 17 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Echogenic tip

ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush




The ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush is used for laparoscopic or ultrasound - guided transvaginal retrieval of oocytes from ovarian follicles.

- Aspiration and flushing line is made possible through a single channel
- Echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With luer lock connector
- Dismantling of luer lock connector
- Ergonomic handle
- 100 cm aspiration line




| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|------------------|--|--------------|--------------------|------------|
| A-OVM-SL16/35-EF | ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush | 16 | 35 | 10 |
| A-OVM-SL17/35-EF | ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush | 17 | 35 | 10 |
| A-OVM-SL18/35-EF | ACE-EF Single Lumen Ovum Pickup Needle with Easy Flush | 18 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Dismantling of
luer lock connector



Echogenic tip

ACE-M Manual Single Lumen Ovum Pickup Needle



The ACE-M Manual Single Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from follicles.

- When used with ultrasonic equipment, echo tip enhances visualization
- Handle designed for greater control during ovum retrieval
- Luer lock connector for connecting a syringe
- The amount of suction applied is manually controlled


Tri-facet geometry
with a echogenic tip



Ergonomic
handle


| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|-----------------|--|--------------|--------------------|------------|
| A-OVM-SL16/35-M | ACE-M Manual Single Lumen Ovum Pickup Needle | 16 | 35 | 10 |
| A-OVM-SL17/35-M | ACE-M Manual Single Lumen Ovum Pickup Needle | 17 | 35 | 10 |
| A-OVM-SL18/35-M | ACE-M Manual Single Lumen Ovum Pickup Needle | 18 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Echogenic tip

DUALFLOW Double Lumen Ovum Pickup Needle



The DUALFLOW OVUMPICK Double Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- The outer lumen allows flushing of follicles




| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|-------------|--|--------------|--------------------|------------|
| DUF-DL17/35 | DUALFLOW Double Lumen Ovum Pickup Needle | 17 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Echogenic tip

OVUMPICK Double Lumen Ovum Pickup Needle



The OVUMPICK Double Lumen Ovum Pickup Needle is used for laparoscopic or ultrasound-guided transvaginal retrieval of oocytes from ovarian follicles.

- Highly echogenic tip
- Silicone bung for a tight fit
- Non-coring needle
- With Luer lock connector
- The outer lumen allows flushing of follicles




| Order No. | Description | Needle Gauge | Needle Length (cm) | Qty. / Box |
|-------------|--|--------------|--------------------|------------|
| OVM-DL16/35 | OVUMPICK Double Lumen Ovum Pickup Needle | 16 | 35 | 10 |
| OVM-DL17/35 | OVUMPICK Double Lumen Ovum Pickup Needle | 17 | 35 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Echogenic tip

EMBRYO TRANS Embryo Transfer Catheter



The EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter with smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately




| Order No. | Description | Catheter Length (cm) | Qty. / Box |
|-----------|---|----------------------|------------|
| EB-ETC23 | EMBRYO TRANS Embryo Transfer Catheter | 23 | 10 |
| EB-ETC18 | EMBRYO TRANS Embryo Transfer Catheter | 18 | 10 |
| STY23 | Stylet for Embryo Transfer Catheter Embryo Trans / Echo Trans | 23 | 10 |
| STY18 | Stylet for Embryo Transfer Catheter Embryo Trans / Echo Trans | 18 | 10 |
| STY23-OS | Stylet for Introducing Catheter 23cm with Outer Sheath | 23 | 10 |
| STY18-OS | Stylet for Introducing Catheter 18cm with Outer Sheath | 18 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Soft inner catheter

ECHO TRANS Embryo Transfer Catheter




The ECHO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic inner catheter facilitates enhanced visibility under ultrasound
- A soft inner catheter with a smooth tip
- Marking on inner catheter and outer sheath
- Shapeable outer sheath
- Stylet for outer sheath available separately




| Order No. | Description | Catheter Length (cm) | Qty. / Box |
|-----------|---|----------------------|------------|
| ET-ETC23 | ECHO TRANS Embryo Transfer Catheter | 23 | 10 |
| ET-ETC18 | ECHO TRANS Embryo Transfer Catheter | 18 | 10 |
| STY23 | Stylet for Embryo Transfer Catheter Embryo Trans / Echo Trans | 23 | 10 |
| STY18 | Stylet for Embryo Transfer Catheter Embryo Trans / Echo Trans | 18 | 10 |
| STY23-OS | Stylet for Introducing Catheter 23cm with Outer Sheath | 23 | 10 |
| STY18-OS | Stylet for Introducing Catheter 18cm with Outer Sheath | 18 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Full echogenic
inner catheter

Trial Transfer Catheter



The Trial Transfer catheter is a closed-ended catheter that is only to be used to determine the passage through the cervix prior to embryo transfer accessible to Allwin's Embryo Transfer Catheter.

- Allows for confident positioning of the outer sheath before introducing the embryo loaded catheter
- Same soft material and smooth tip as catheter to minimize patient trauma
- Centimetre depth marking for confident and precise catheter placement
- Limulus Amebocyte Lysate (LAL) tested to ensure the product is non-pyrogenic




| Order No. | Description | Catheter Length (cm) | Qty. / Box |
|-------------|---|----------------------|------------|
| TT-EB-ETC23 | Trial Transfer Catheter for Embryo Trans 23cm | 23 | 10 |
| TT-EB-ETC18 | Trial Transfer Catheter for Embryo Trans 18cm | 18 | 10 |
| TT-ET-ETC23 | Trial Transfer Catheter for Echo Trans 23cm | 23 | 10 |
| TT-ET-ETC18 | Trial Transfer Catheter for Echo Trans 18cm | 18 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip




The BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe and outer sheath
- Bulb shaped tip of the outer sheath
- Silicone depth positioner for precise placement
- The bulb-shaped angled outer sheath eases passage through the cervix




| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| BT-ETC-U | BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip | 10 |
| STY-BU | Stylet for BULB TRANS ULTRA Embryo Transfer Catheter - Bulb Tip | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Bulb tip with
soft inner catheter

BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip





The BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic band on the inner catheter to provide visibility under ultrasound
- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe and outer sheath
- Bulb shaped angled tip of the outer sheath
- Silicone depth positioner for precise placement
- The catheter is curved to facilitate insertion



| Order No. | Description | Qty. / Box |
|------------|--|------------|
| BT-ETC-U-S | BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip | 10 |
| STY-BU | Stylet for BULB TRANS STAR Embryo Transfer Catheter - Bulb Tip | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Inner catheter tip
echogenic metal band

BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip



The BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe
- Bulb shaped tip of the outer sheath
- Silicone depth positioner for precise placement
- The bulb-shaped angled outer sheath eases passage through the cervix




| Order No. | Description | Qty. / Box |
|------------|--|------------|
| BT-ETC-U-P | BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip | 10 |
| STY-BU | Stylet for BULB TRANS ULTRA PLUS Embryo Transfer Catheter - Bulb Tip | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Premium bulb tip catheter
with soft inner catheter

BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip





The BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Firm outer sheath
- A supported soft inner catheter with a smooth tip
- Marking on inner steel pipe
- Bulb shaped angled tip of the outer sheath
- Silicone depth positioner for precise placement
- The bulb-shaped angled outer sheath eases passage through the cervix
- Echogenic band on the inner catheter to improve the visibility under ultrasound



| Order No. | Description | Qty. / Box |
|--------------|---|------------|
| BT-ETC-U-S-P | BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip | 10 |
| STY-BU | Stylet for BULB TRANS STAR PLUS Embryo Transfer Catheter - Bulb Tip | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Inner catheter tip
echogenic metal band

PRO EMBRYO TRANS Embryo Transfer Catheter



The PRO EMBRYO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Silicone depth positioner for precise placement
- Orientation arrow on the hub
- Angled outer sheath provides easy passage through the cervical canal
- Available in compatible stylet or soft obturator



Soft Obturator for PRO EMBRYO TRANS Embryo Transfer Catheter



Stylet for PRO EMBRYO TRANS Embryo Transfer Catheter


| Order No. | Description | Qty./Box |
|-------------|---|----------|
| PEB-ETC | PRO EMBRYO TRANS Embryo Transfer Catheter | 10 |
| PEB-ETC-SB | PRO EMBRYO TRANS Embryo Transfer Catheter with Soft Obturator | 10 |
| PEB-ETC-STY | PRO EMBRYO TRANS Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Stylet for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Supported inner catheter
with smooth tip

PRO EMBRYO STAR Embryo Transfer Catheter



The PRO EMBRYO STAR Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic band on the inner catheter to provide visibility under ultrasound
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Silicone depth positioner for precise placement
- Orientation arrow on the hub
- Angled outer sheath provides easy passage through the cervical canal
- Available in compatible stylet or soft obturator




Soft Obturator for PRO EMBRYO STAR Embryo Transfer Catheter




Stylet for PRO EMBRYO STAR Embryo Transfer Catheter

| Order No. | Description | Qty./Box |
|---------------|---|----------|
| PEB-ETC-S | PRO EMBRYO STAR Embryo Transfer Catheter | 10 |
| PEB-ETC-S-SB | PRO EMBRYO STAR Embryo Transfer Catheter with Soft Obturator | 10 |
| PEB-ETC-S-STY | PRO EMBRYO STAR Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Stylet for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter | 10 |



Echogenic metal band

PRO ECHO TRANS Embryo Transfer Catheter



The PRO ECHO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- Echogenic inner catheter facilitates enhanced visibility under ultrasound
- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Silicone depth positioner for precise placement
- Orientation arrow on hub
- Angled outer sheath provides easy passage through the cervical canal
- Available in compatible stylet or soft obturator




Soft Obturator for PRO ECHO TRANS Embryo Transfer Catheter



Stylet for PRO ECHO TRANS Embryo Transfer Catheter


| Order No. | Description | Qty. / Box |
|-------------|---|------------|
| PET-ETC | PRO ECHO TRANS Embryo Transfer Catheter | 10 |
| PET-ETC-SB | PRO ECHO TRANS Embryo Transfer Catheter with Soft Obturator | 10 |
| PET-ETC-STY | PRO ECHO TRANS Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Stylet for Pro Embryo Trans / Pro Embryo Star / Pro Echo Trans Embryo Transfer Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Full echogenic
inner catheter

PRO PASS Embryo Transfer Catheter



The PRO PASS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Orientation arrow on the hub
- Available in compatible stylet or soft obturator

Soft inner catheter
and smooth tip

Centimetre marking

Orientation arrow on hub



Soft Obturator for PRO PASS Embryo Transfer Catheter



Stylet for PRO PASS Embryo Transfer Catheter




| Order No. | Description | Qty. / Box |
|------------|--|------------|
| PP-ETC | PRO PASS Embryo Transfer Catheter | 10 |
| PP-ETC-SB | PRO PASS Embryo Transfer Catheter with Soft Obturator | 10 |
| PP-ETC-STY | PRO PASS Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Stylet for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



PRO PASS STAR Embryo Transfer Catheter



The PRO PASS STAR Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Orientation arrow on the hub
- Available in compatible stylet or soft obturator





Soft Obturator for PRO PASS STAR Embryo Transfer Catheter



Stylet for PRO PASS STAR Embryo Transfer Catheter


| Order No. | Description | Qty. / Box |
|--------------|--|------------|
| PP-ETC-S | PRO PASS STAR Embryo Transfer Catheter | 10 |
| PP-ETC-S-SB | PRO PASS STAR Embryo Transfer Catheter with Soft Obturator | 10 |
| PP-ETC-S-STY | PRO PASS STAR Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Stylet for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Full echogenic inner catheter

PRO PASS ECHO TRANS Embryo Transfer Catheter



The PRO PASS ECHO TRANS Embryo Transfer Catheter is used to introduce in-vitro fertilized (IVF) embryos into the uterine cavity.

- A soft inner catheter & smooth tip
- Marking on inner steel pipe & outer sheath
- Flexible catheter
- Orientation arrow on the hub
- Available in compatible stylet or soft obturator

Visible under
ultrasound

Centimetre marking

Orientation arrow on hub




Soft Obturator for PRO PASS ECHO TRANS Embryo Transfer Catheter


Stylet for PRO PASS ECHO TRANS Embryo Transfer Catheter

| Order No. | Description | Qty. / Box |
|-------------|--|------------|
| PPE-ETC | PRO PASS ECHO TRANS Embryo Transfer Catheter | 10 |
| PPE-ETC-SB | PRO PASS ECHO TRANS Embryo Transfer Catheter with Soft Obturator | 10 |
| PPE-ETC-STY | PRO PASS ECHO TRANS Embryo Transfer Catheter with Stylet | 10 |
| SB-PE | Soft Obturator for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Stylet for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Trial Transfer Catheter for Pro Pass



The Trial Transfer Catheter for Pro Pass is a closed-ended catheter that is only to be used to determine the passage through the cervix prior to embryo transfer accessible to Allwin's Embryo Transfer Catheter.

- Allows for confident positioning of the outer sheath before introducing the embryo loaded catheter
- Same soft material and smooth tip as catheter to minimize patient trauma
- Centimetre depth marking for confident and precise catheter placement
- Limulus Amebocyte Lysate (LAL) tested to ensure the product is non-pyrogenic




| Order No. | Description | Qty. / Box |
|------------|--|------------|
| TT-PP-ETC | Trial Transfer Catheter for Pro Pass Embryo Transfer Catheter | 10 |
| TT-PPS-ETC | Trial Transfer Catheter for Pro Pass Star Embryo Transfer Catheter | 10 |
| TT-PPE-ETC | Trial Transfer Catheter for Pro Pass Echo Embryo Transfer Catheter | 10 |
| SB-PE | Soft Obturator for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |
| STY-PE | Styler for Pro Pass / Pro Pass Star / Pro Pass Echo Trans Embryo Transfer Catheter | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




INSEMINA Intra Uterine Insemination Catheter



The INSEMINA Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Open tip as well as closed tip (with two lateral holes) catheters available



| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| I-IUI | INSEMINA Intra Uterine Insemination Catheter | 25 |
| I-IUI-O | INSEMINA Intra Uterine Insemination Catheter (Open Tip) | 25 |

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip




Open Tip

CURVED Intra Uterine Insemination Catheter



The CURVED Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Curve shaped catheter
- Open tip as well as closed tip (with two lateral holes) catheters available



| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| I-IUI-C | CURVED Intra Uterine Insemination Catheter | 25 |
| I-IUI-C-O | CURVED Intra Uterine Insemination Catheter (Open Tip) | 25 |

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip




Open Tip

TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet



The TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- A soft uterine catheter with firm support
- Integrally molded hub
- Silicone depth positioner for precise placement
- Malleable stylet to aid catheter insertion
- Open tip as well as closed tip (with two lateral holes) catheters available



| Order No. | Description | Qty. / Box |
|-----------|--|------------|
| T-IUI | TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet | 25 |
| T-IUI-O | TRACKABLE Intra Uterine Insemination Catheter with Malleable Stylet (Open Tip) | 25 |

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip



Open Tip

FLEXI Intra Uterine Insemination Catheter



The FLEXI Intra Uterine Insemination Catheter (Open Tip & Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- Shapeable outer sheath for easy insertion
- Centimetre marking on a catheter for precise placement
- Stylet available separately
- Open tip as well as closed tip (with two lateral holes) catheters available

Centimetre marking provided



Soft uterine catheter

Shapeable
outer sheath

Firm
support



Stylet

| Order No. | Description | Qty. / Box |
|-----------|--|------------|
| F-IUI | FLEXI Intra Uterine Insemination Catheter | 25 |
| F-IUI-O | FLEXI Intra Uterine Insemination Catheter (Open Tip) | 25 |
| STY18 | Stylet for FLEXI IUI Intra Uterine Insemination Catheter | 25 |

*All products may not be available in all countries. Please contact our sales department for details.




Close Tip




Open Tip

INTEGRA Intra Uterine Insemination Catheter



The INTEGRA Intra Uterine Insemination Catheter (Close Tip) is intended for the introduction of washed spermatozoa into the uterine cavity.

- Placement becomes easier with the distal tip's soft, flexible curve.
- The stainless steel inner cannula is malleable.
- Luer lock connector.
- Silicone stopper for adjustable positioning.



| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| IG-IUI | INTEGRA Intra Uterine Insemination Catheter | 25 |

*All products may not be available in all countries. Please contact our sales department for details.



Side Hole

STUMPY Intra Uterine Insemination Cannula



The STUMPY Intra Uterine Insemination Cannula is intended for the introduction of washed spermatozoa into the uterine cavity.

- Rounded tip
- Luer lock
- Flexible cannula
- Small size catheter




| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| ST-IUI | STUMPY Intra Uterine Insemination Cannula | 25 |

*All products may not be available in all countries. Please contact our sales department for details.




SALPINX CATH Fallopian Tube Recanalization Catheter



The SALPINX CATH Fallopian Tube Recanalization Catheter is used for fallopian tube recanalization under radiological fluoroscopic control.



- Preshaped 5 Fr catheter
- Centimetre marking on the catheter for precise placement
- Malleable stylet
- Hydrophilic guidewire available separately



| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| SC-FTR-M | SALPINX CATH Fallopian Tube Recanalization Catheter with Malleable Stylet | 5 |
| SC-FTR-MG | SALPINX CATH Fallopian Tube Recanalization Catheter with Malleable Stylet and Hydrophilic Guidewire | 1 |


*All products may not be available in all countries. Please contact our sales department for details.

HYSTERO CATH Hysterosalpingography Catheter



The HYSTERO CATH Hysterosalpingography Catheter is used for injection of contrast media into the uterine cavity and fallopian tube(s) to evaluate the patency of the fallopian tubes.

- Strong low profile 7 Fr balloon
- Closed tip with side hole and open tip HSG catheter available
- Stopcock designed for locking the balloon channel and regulating the flow of contrast media.
- Stylet provided for added rigidity during the introduction
- 3 ml syringe included





| Order No. | Description | Qty. / Box |
|-----------|--|------------|
| HC-HSG7 | HYSTERO CATH Hysterosalpingography Catheter | 5 |
| HC-HSG-O7 | HYSTERO CATH Hysterosalpingography Catheter Open Taper Tip | 5 |

*All products may not be available in all countries. Please contact our sales department for details.

Caution:


This product contains natural rubber latex which may cause allergic reaction in some individuals.

HYSTERO CATH SILICONE Hysterosalpingography Catheter



The HYSTERO CATH SILICONE Hysterosalpingography Catheter is used for injection of contrast media into the uterine cavity and fallopian tube(s) to evaluate the patency of the fallopian tubes.


- Strong low profile 7 Fr balloon
- Closed tip with side hole and open tip HSG catheter available
- Stopcock designed for locking the balloon channel and regulating the flow of contrast media.
- Stylet provided for added rigidity during the introduction
- 3 ml syringe included



| Order No. | Description | Qty. / Box |
|--------------|---|------------|
| HCSIL-HSG-7 | HYSTERO CATH SILICONE Hysterosalpingography Catheter | 5 |
| HCSIL-HSG-O7 | HYSTERO CATH SILICONE Hysterosalpingography Catheter Open Taper Tip | 5 |

*All products may not be available in all countries. Please contact our sales department for details.

THE SAMPLER Endometrial Biopsy Catheter



THE SAMPLER Endometrial Biopsy Catheter is used for performing histological tissue biopsy and extracting samples of uterine menstrual content. Used for obtaining endometrial tissue.

- Centimetre marking provided on the catheter for precise placement
- Smooth rounded tip

Smooth rounded tip




Centimetre marking

| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| SEBC | THE SAMPLER Endometrial Biopsy Catheter | 25 |

*All products may not be available in all countries. Please contact our sales department for details.

PRO-VAC Endometrial Biopsy Sampler



The PRO-VAC Endometrial Biopsy Sampler is used for performing histological tissue biopsy and extracting samples of uterine menstrual content. Used for obtaining endometrial tissue.

- Centimetre marking provided on the catheter for precise placement
- Smooth rounded curved tip
- Locking syringe

Smooth rounded curved tip




Centimeter marking provided

Locking syringe

| Order No. | Description | Qty. / Box |
|-----------|------------------------------------|------------|
| PVES | PRO-VAC Endometrial Biopsy Sampler | 25 |

*All products may not be available in all countries. Please contact our sales department for details.

SURESTRIP Vitrification Straws



The SURESTRIP Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- 1 - cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours

Ring marking for easy visualization during insertion

Transparent cap



Marking on tip for easy visibility under microscope

Flat surface wide body for writing


Strip width ensuring minimum placement of solution

| Colour | White | Red | Ivory | Blue | Yellow | Qty. / Box |
|-----------|-------|-------|-------|-------|--------|------------|
| Order No. | SUS-W | SUS-R | SUS-I | SUS-B | SUS-Y | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



SURESTRIP Vitrification Straws - Weighted



The SURESTRIP Vitrification Straws-Weighted is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- Weighted balls prevents it from floating in nitrogen
- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- 1 - cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours

Ring marking for easy visualization during insertion

Transparent cap

Weighted balls



Marking on tip for easy visibility under microscope



Flat surface wide body for writing

Strip width ensuring minimum placement of solution


| Colour | White | Red | Ivory | Blue | Yellow | Qty. / Box |
|-----------|----------|----------|----------|----------|----------|------------|
| Order No. | SUS-W-WT | SUS-R-WT | SUS-I-WT | SUS-B-WT | SUS-Y-WT | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



Weighted balls

SURESTRIP Vitrification Straws - Weighted Plug



The SURESTRIP Vitrification Straws - Weighted Plug is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- Weight prevents it from floating in nitrogen
- 1 - cell MEA $\geq 80\%$ expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours
- One goblet is filled with eight vitrification straws

Ring marking for easy visualization during insertion

Transparent cap

Weight

Plug



Marking on tip for easy visibility under microscope

Flat surface wide body for writing

Strip width ensuring minimum placement of solution




| Colour | White | Red | Ivory | Blue | Yellow | Qty. / Box |
|-----------|-----------|-----------|-----------|-----------|-----------|------------|
| Order No. | SUS-W-WTP | SUS-R-WTP | SUS-I-WTP | SUS-B-WTP | SUS-Y-WTP | 10 |

*All products may not be available in all countries. Please contact our sales department for details.




Weighted Plug

SURELOCK Vitrification Straws



The SURELOCK Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- Concave tip for exact and secure placement of specimen under microscope
- Secure air tight coupling
- Semi-transparent body and cap
- The arrow on top surface ensuring orientation during use
- Easy visibility under microscope
- 1 - cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of 1-cell stage of embryos
- Available in 5 colours




| Colour | Pink | Yellow | Green | Blue | Violet | Qty. / Box |
|-----------|-------|--------|-------|-------|--------|------------|
| Order No. | SUL-P | SUL-Y | SUL-G | SUL-B | SUL-V | 10 |

*All products may not be available in all countries. Please contact our sales department for details.



D-TIP Denudation Pipettes






The Device is intended to be used for denudation i.e. to remove cumulus oocytes complexes from oocytes before procedure, "Intracytoplasmic Sperm Injection (ICSI)" and handling or manipulating of oocytes and embryos in Assisted reproduction technique (ART).

- Designed to adapt to surfaces with flexibility, lowering the risk of scratching or breaking
- Precisely engineered for delicate and accurate handling
- Different Color code for different sizes of pipette
- Endotoxin Tested (LAL test) 20EU/device
- Mouse Embryo Assay (MEA) $\geq 80\%$
- Gamma Sterilized



| Order No. | Description | Colour | Qty/ Vial | Vial / Box |
|-------------|---------------------------------|--------|-----------|------------|
| D-TIP 135µm | D-TIP Denudation Pipettes 135µm | Yellow | 10 | 5 |
| D-TIP 145µm | D-TIP Denudation Pipettes 145µm | White | 10 | 5 |
| D-TIP 155µm | D-TIP Denudation Pipettes 155µm | Purple | 10 | 5 |
| D-TIP 175µm | D-TIP Denudation Pipettes 175µm | Red | 10 | 5 |
| D-TIP 275µm | D-TIP Denudation Pipettes 275µm | Gray | 10 | 5 |
| D-TIP 300µm | D-TIP Denudation Pipettes 300µm | Brown | 10 | 5 |

*All products may not be available in all countries. Please contact our sales department for details.



AspiPump 90D

- High vacuum, low flow pump capable of operating at normal vacuum levels ranging from -10mmHg to -300mmHg and boost vacuum levels from -300mmHg to -500mmHg, with adjustable vacuum settings in increments of -1mmHg.
- Red LED display for easy readability.
- Equipped with a dual foot pedal featuring color coding: Yellow for boost operation and blue for normal operation.
- The Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors (200cm), and Aspi Hydrophobic Filter Set are provided sterile & intended for single use only.

Aspi Suction Bottle with Connector



**COMING
SOON**

ThermiSure TBT-12
(Warmer)

Aspi Hydrophobic Filter Set



Aspi Suction Line with
Luer Lock Connectors




AspiPump Foot Switch



| Order No. | Description | Qty. / Box |
|-----------|---|------------|
| ASPI 90D | AspiPump 90D Low flow / high vacuum DC operated pump, dual foot pedal, Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors & Aspi Hydrophobic Filter Set, Universal adaptor, User guide Manual. | 1 |
| ASP-SB | Aspi Suction Bottle with Connector Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units. | 4 |
| ASP-HFS | Aspi Hydrophobic Filter Set Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units. | 4 |
| ASP-SL | Aspi Suction Line with Luer Lock Connectors Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units. | 4 |
| ASP-FS | AspiPump Foot Switch Suitable for Allwin AspiPump 90D. | 4 |


*All products may not be available in all countries. Please contact our sales department for details.

ASPI Aspiration Needle



The ASPI Aspiration Needle is used for laparoscopic or ultrasound-guided transvaginal aspiration.


- Bevel tip
- Centimetre marking provided for precise placement
- Transparent luer lock connector



| Order No. | Description | Gauge | Size cm | Qty. / Box |
|------------|------------------------|-------|---------|------------|
| ASPI 18/30 | ASPI Aspiration Needle | 18 | 30 | 10 |
| ASPI 22/30 | ASPI Aspiration Needle | 22 | 30 | 10 |

*All products may not be available in all countries. Please contact our sales department for details.

Female Reproductive System





Medical Devices

With Allwin, We All Win



US Sales

✉ sales.usa@allwinmedical.com

☎ +1 714-572-1709

India Sales

✉ sales.india@allwinmedical.com

☎ +91 99870-84563

Worldwide Sales (Excluding India/US)

✉ sales@allwinmedical.com

☎ +91 95949-81581

🌐 www.allwinmedical.com