


SURESTRIP Vitrification Straws



The SURESTRIP Vitrification Straw is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- 1 - cell MEA $\geq 80\%$ expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours

Ring marking for easy visualization during insertion

Transparent cap



Marking on tip for easy visibility under microscope

Flat surface wide body for writing

Strip width ensuring minimum placement of solution



Colour	Pink	Red	Green	Blue	Yellow	Qty. / Box
Order No.	SUS-P	SUS-R	SUS-G	SUS-B	SUS-Y	10

*All products may not be available in all countries. Please contact our sales department for details.

