



Glider Amplatz Sheath



The Glider Amplatz Sheath is used to maintain previously established nephrostomy tract.

- Smooth surface minimizes trauma to patient
- Bevel tip designed for atraumatic placement and ease of advancement
- Firm and kink-resistant
- Radiopaque for easier visualization
- Available in 12 Fr to 30 Fr even sizes
- Smooth surface assists in easier introduction of catheter



Order No.	Fr	Length (cm)	Qty. / Box
AM12	12	17	5
AM14	14	17	5
AM16	16	17	5
AM18	18	17	5
AM20	20	17	5
AM22	22	17	5
AM24	24	17	5
AM26	26	17	5
AM28	28	17	5
AM30	30	17	5

*All products may not be available in all countries. Please contact our sales department for details.