



## D-TIP Denudation Pipettes



The Device is intended to be used for denudation i.e. to remove cumulus oocytes complexes from oocytes before procedure, "Intracytoplasmic Sperm Injection (ICSI)" and handling or manipulating of oocytes and embryos in Assisted reproduction technique (ART).

- Designed to adapt to surfaces with flexibility, lowering the risk of scratching or breaking
- Precisely engineered for delicate and accurate handling
- Different Color code for different sizes of pipette
- Endotoxin Tested (LAL test) 20EU/device
- Mouse Embryo Assay (MEA)  $\geq 80\%$
- Gamma Sterilized



| Order No.   | Description                     | Colour | Qty/ Vial | Vial / Box |
|-------------|---------------------------------|--------|-----------|------------|
| D-TIP 135µm | D-TIP Denudation Pipettes 135µm | Yellow | 10        | 5          |
| D-TIP 145µm | D-TIP Denudation Pipettes 145µm | White  | 10        | 5          |
| D-TIP 155µm | D-TIP Denudation Pipettes 155µm | Purple | 10        | 5          |
| D-TIP 175µm | D-TIP Denudation Pipettes 175µm | Red    | 10        | 5          |
| D-TIP 275µm | D-TIP Denudation Pipettes 275µm | Grey   | 10        | 5          |
| D-TIP 300µm | D-TIP Denudation Pipettes 300µm | Brown  | 10        | 5          |

\*All products may not be available in all countries. Please contact our sales department for details.

