

ASPI Aspiration Needle

Description:

ASPI Aspiration Needle

Intended use:

Aspiration needle is used for aspiration of laparoscopic or ultrasound guided transvaginal aspiration.

Contraindications:

There are no known contraindications.

Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician.

Warning / Precaution:

- Do not use if package is opened or damaged.
- Only experienced surgeon should place this device.
- This device is supplied sterile and intended for single use only.
- Attempt to resterilize a single use device may compromise its structural integrity and/or lead to device failure which may result in patient injury or illness.
- Storage temperature between 5°C to 30°C.
- Dispose in accordance with recognized medical practice and under observance of applicable law and regulation.

Instruction for Use:

Using fluoroscopy / ultrasound guidance, advance the needle through anesthetized cut fascia, till it reaches the point of interest.

Further steps should be carried out by the physician as per requirement.

Do not reuse, reprocess & resterilize:

Reuse, Reprocessing & resterilization may compromise the structural integrity, can also create risk of contamination & may not give desired result or create complications, infections which may result in injury, illness or death.

Limited Express Warranty:

The limited express warranty as set forth herein is exclusive and in lieu of all warranties of merchantability and fitness for use, remedies, obligations and liabilities for consequential damages. The products are being sold only for the purpose described herein and such limited express warranty runs only to the original user. In no event shall ALLWIN be liable for any breach of warranty in any amount exceeding the purchase price of the product. ALLWIN reserves the right to make design changes to products without liability to incorporate said changes in ALLWIN products previously designed or sold.

