

## AspiPump 90D

- High vacuum, low flow pump capable of operating at normal vacuum levels ranging from -10mmHg to -300mmHg and boost vacuum levels from -300mmHg to -500mmHg, with adjustable vacuum settings in increments of -1mmHg.
- Red LED display for easy readability.
- Equipped with a dual foot pedal featuring color coding: Yellow for boost operation and blue for normal operation.
- The Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors (200cm), and Aspi Hydrophobic Filter Set.

Aspi Suction Bottle with Connector



ThermiSure T12



Aspi Hydrophobic Filter Set



Aspi Suction Line with Luer Lock Connectors



Aspi Foot Switch



Order No.	Description	Qty. / Box
ASPI 90D	<b>AspiPump 90D</b> Low flow / high vacuum DC operated pump, dual foot pedal, Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors & Aspi Hydrophobic Filter Set, Universal adaptor.	1
ASP-SB	<b>Aspi Suction Bottle with Connector</b> Supplied with Allwin AspiPump 90D, Provided in cartons of 4 units.	4
ASP-HFS	<b>Aspi Hydrophobic Filter Set</b> Supplied with Allwin AspiPump 90D, Provided in cartons of 4 units.	4
ASP-SL	<b>Aspi Suction Line with Luer Lock Connectors</b> Supplied with Allwin AspiPump 90D, Provided in cartons of 4 units.	4
ASP-FS	<b>Aspi Foot Switch</b> Suitable for Allwin AspiPump 90D.	1

\*All products may not be available in all countries. Please contact our sales department for details.