



## D-TIP Denudation Pipettes



The Device is intended to be used for denudation i.e. to remove cumulus oocytes complexes from oocytes before procedure, "Intracytoplasmic Sperm Injection (ICSI)" and handling or manipulating of oocytes and embryos in Assisted reproduction technique (ART).

- Designed to adapt to surfaces with flexibility, lowering the risk of scratching or breaking
- Precisely engineered for delicate and accurate handling
- Different Color code for different sizes of pipette
- Endotoxin Tested (LAL test) 20EU/device
- Mouse Embryo Assay (MEA)  $\geq 80\%$
- Gamma Sterilized



Order No.	Description	Colour	Qty/ Vial	Vial / Box
D-TIP 125µm	<b>D-TIP</b> Denudation Pipettes 125µm	Blue	10	5
D-TIP 130µm	<b>D-TIP</b> Denudation Pipettes 130µm	Yellow	10	5
D-TIP 135µm	<b>D-TIP</b> Denudation Pipettes 135µm	Yellow	10	5
D-TIP 140µm	<b>D-TIP</b> Denudation Pipettes 140µm	White	10	5
D-TIP 145µm	<b>D-TIP</b> Denudation Pipettes 145µm	White	10	5
D-TIP 155µm	<b>D-TIP</b> Denudation Pipettes 155µm	Purple	10	5
D-TIP 170µm	<b>D-TIP</b> Denudation Pipettes 170µm	Red	10	5
D-TIP 175µm	<b>D-TIP</b> Denudation Pipettes 175µm	Red	10	5
D-TIP 200µm	<b>D-TIP</b> Denudation Pipettes 200µm	Pink	10	5
D-TIP 275µm	<b>D-TIP</b> Denudation Pipettes 275µm	Grey	10	5
D-TIP 300µm	<b>D-TIP</b> Denudation Pipettes 300µm	Brown	10	5

\*All products may not be available in all countries. Please contact our sales department for details.

