

## AspiPump 90D

- High vacuum, low flow pump capable of operating at normal vacuum levels ranging from -10mmHg to -300mmHg and boost vacuum levels from -300mmHg to -500mmHg, with adjustable vacuum settings in increments of -1mmHg.
- Red LED display for easy readability.
- Equipped with a dual foot pedal featuring color coding: Yellow for boost operation and blue for normal operation.
- The Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors (200cm), and Aspi Hydrophobic Filter Set are provided sterile & intended for single use only.

Aspi Suction Bottle with Connector



ThermiSure T12



Aspi Hydrophobic Filter Set



Aspi Suction Line with Luer Lock Connectors



Aspi Foot Switch



| Order No. | Description                                                                                                                                                                                                                                   | Qty. / Box |
|-----------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| ASPI 90D  | <b>AspiPump 90D</b><br>Low flow / high vacuum DC operated pump, dual foot pedal, Aspi Suction Bottle with Connector (100ml), Aspi Suction Line with Luer Lock Connectors & Aspi Hydrophobic Filter Set, Universal adaptor, User guide Manual. | 1          |
| ASP-SB    | <b>Aspi Suction Bottle with Connector</b><br>Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units.                                                                                              | 4          |
| ASP-HFS   | <b>Aspi Hydrophobic Filter Set</b><br>Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units.                                                                                                     | 4          |
| ASP-SL    | <b>Aspi Suction Line with Luer Lock Connectors</b><br>Supplied with Allwin AspiPump 90D, Provided sterile, for single use only and in cartons of 4 units.                                                                                     | 4          |
| ASP-FS    | <b>Aspi Foot Switch</b><br>Suitable for Allwin AspiPump 90D.                                                                                                                                                                                  | 1          |

\*All products may not be available in all countries. Please contact our sales department for details.