


SURESTRIP Vitrification Straws - Weighted Plug



The SURESTRIP Vitrification Straws - Weighted Plug is a cryopreservation storage device. It is used in the vitrification procedure to contain and maintain oocytes and embryos.

- High survival and development rates
- Secure coupling
- Marking on tip extreme for easy visibility under the microscope
- Adequate flat surface for writing
- Tip width ensuring placement of minimum solution
- Weight prevents it from floating in nitrogen
- 1 - cell MEA $\geq$ 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma sterilized
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours
- One goblet is filled with eight vitrification straws



Colour	Pink	Red	Green	Blue	Yellow	Qty. / Box
Order No.	SUS-P-WTP	SUS-R-WTP	SUS-G-WTP	SUS-B-WTP	SUS-Y-WTP	10

*All products may not be available in all countries. Please contact our sales department for details.



Weighted Plug